

# Single Use Cranial Access Kit



Cost effective and efficient. Each kit is supplied as a sterile pack

Kit Type

Order Code

UOM

**Standard Kit**  
(detailed overleaf)

**RMT0082**

**10**

**Standard Kit + CH5 Drill Burr**

**RMT0082-CH5**

**10**

Cranial Hand Drill

Order Code


UOM

**Plastic Cranial Hand Drill**

**D-5001-S**

**1**

Drills are compatible with cranial access burrs from all manufacturers



**DELTA**

Delta Surgical Ltd. ☎ 01782 637009 ✉ sales@deltasurgical.co.uk  
Unit 10, Evolution, Lymedale Business Park, Newcastle-under-Lyme,  
Staffordshire, ST5 9QF



# Standard Kit Contents

This convenient, pre-packaged sterile kit provides access to the subarachnoid space. Designed to be used with external drainage and monitoring to reduce intracranial pressure. Each kit contains the following:

- Holzheimer Blunt 10cm
- Disposable Hand Drill
- Mosquito Forceps, Curved 13cm
- Mayo Needle Holder, 16cm
- Adson Tissue Forceps 12cm
- Adson Dressing Forceps
- Iris Scissor, 11cm
- Bone Wax
- Suture Mersilk 3/8 circle RC 26mm, black 2/0, 45cms
- Suture Ethilon W1625
- Skin Marker Pen, Fine (with Ruler)
- Syringe, 10ml LS Ecc Emerald
- Bard Parker Handle with No 11 Blade
- Bard Parker Handle with No 15 Blade
- Sponge Stick 20cm Orange
- Hand towel ECO 45x33cm
- Swab Gauze XRD 7.5 x 7.5cm 8ply
- Tray Full 250x290x50mm
- Tray Half Shallow 245x135x20mm
- Drape 75 x 75cm Fen 6 x 8cm S/A SMS

A selection of the items listed above are displayed below. For more details or to place an order please contact us: [sales@deltasurgical.co.uk](mailto:sales@deltasurgical.co.uk)  
01782 637009

